

BioModeling Sand (B-MDL1SD)

Check: Well shakened resin, Clean optical pathways, Good resin tank condition.

Recommendation:

Prevention of cross-contamination with other non-bio resins ensures optimal biocompatibility for printed models. Dedicate accessories for BioModeling Sand: resin tank, build platform and wash.

Print Settings:

| Brand | Printer Model | Settings |
|---|--|--|
| Open system MSLA <ul style="list-style-type: none"> ● Elegoo ● Phrozen ● Anycubic | <ul style="list-style-type: none"> ● Mars 5 ultra, Saturn 4...etc. ● Mighty Revo, Sonic mini 8ks, ...etc. ● Photon M3 Premium ...etc. | To get the latest printing settings, simply visit the ApplyLabWork Support page and download the files for your printer model. |
| Asiga | <ul style="list-style-type: none"> ● Pro 4K UV385 ● Max UV385 | |
| HeyGears | Reflex | PARP10 (100um) |

Note: Each printer is unique, and light intensity can vary even among identical models. Minor adjustments may be necessary, especially true for open-system MSLA printers.

Cleaning: [fill fresh 95% IPA in two baths]

A-bath: Submerge and agitate printed models in IPA. Use toothbrush to gently remove uncured resin quickly on the model surfaces. Transfer cleaned model to B-bath.

B-bath: Submerge cleaned models from A-bath for **3 mins**.

Total IPA clean time should not exceed 5 minutes.

Drying:

Blow dry the models with compressed air, let models **sit dry 25 minutes minimum** in a shaded airy place. Printed parts must be free of IPA before UV post-curing. Extended wash time and residual IPA (if not dried completely) causes deformation and lower mechanical properties.

Biocompatible Resin

UV Post-Curing Setting(Required):

| Brand | Post curing Chamber Model | Settings |
|----------|---------------------------|----------------------|
| Elegoo | Mercury Cure | 45 min |
| Formlabs | FormCure | 35°C 30min |
| Formlabs | Fast Cure | Intensity=1, 5min |
| NK Optik | Otoflash | 2500 flashes |
| HeyGears | UltraCraft Cure | PAP10 |

Allow the models to cool naturally to room temperature after curing. (Avoid external forces, no stacking of models during cooling.)

Storage:

Filter out all debris in used-resin before storage. Keep resin in the original bottle, storage temp (15-26°C / 59-79°F) in a dry and dark area.

Note: Pigment settlement is normal. Gently mix remaining resin in resin tank with soft “resin wiper” for color evenness.