


## DLP385 Modeling – Tan (D385M-R001TN) Asiga Max UV385


### Download ini file

<https://www.applylabwork.com/support/>


Step 1. Open Composer software, select **New**.



Step 2. Click  icon, and click **Preferences...**



Step 3. Click  icon, and choose the folder where you download our .ini file, click **OK**.



### Ready to Print:

Choose **ALW-DLP385-Tan** from **Materials** drop-down menu, and choose appropriate layer thickness for your model [high layer thickness reduces printing time but decreases print quality]

Load 3D model and double check support and orientation are optimal before print.

### Washing

Submerge print in clean IPA bath for **5 minutes**, blow off excess IPA, set print in a shaded airy place to dry before curing.

### Curing

**Asiga Flash 10-30 minutes**

Optimal rigidity can be achieved with 30 minutes, reduced curing time will provide more material flexibility.

### Storage and Handling

- Don't leave resin in the tray for more than 3 days without printing.
- Not accessible to Children
- Filter resin in case of fail prints.
- Slight pigment settlement is normal. Gently mix remaining resin in resin tray with **rubber spatula** for color evenness.

Note:

This resin is formulated on **Asiga Max UV385**.